

EZ Ray

Premium Intra-oral X-ray

Dental Pioneer
VATECH
VATECH Global



HealthCare Dental Limited
凱健牙科有限公司
23/F, 8 Commercial Tower, 8 Sun Yip Street, Chai Wan, HKSAR
香港柴灣新業街8號八號商業廣場23樓
Tel: (852) 2563 5193
Fax: (852) 2565 9537, 2805 5931
Email: sales.hcdental@hcdh.com.hk


Designed for your convenience

- Streamline design for easy one-hand grip
- 4.7kg head for stable positioning
- Only 2 clicks needed for capturing
- Detachable control panel




Detachable control panel
Hazard free, Rounded edges

Easy one-hand grip



Product highlights

Safe & Clear Image



EzRay Premium provides diagnostic values with advanced technology.

A focal spot of 0.4 mm create sharpen images perfectly matched to whatever digital intra-oral sensor, PSP or film product you use.

User friendly

► Easy to use

Intuitive icons on LED screen Detachable panel


► Optimized automatic settings

EzRay Premium has a pre-programmed setting for each tooth. With one click of a button, it automatically optimizes image capturing such as capture area and patient condition for each case.




Dimensions

Retracted Wall Mount



Extended Position Wall Mount



Specification (VEX-S100W)

Focal Spot	0.4mm
Tube Voltage	Selectable 50 ~ 70kVp
Current	4 ~ 7 mA
Exposure Time	Selectable 0.04 ~ 2.0 sec
Arm Length (Option)	450, 600, 900 mm
Line Voltage	100 ~ 230V ± 10%
Accessories	Exposure hand switch
Cone	Round 20cm [8"] Option: Round 30cm [12"], Rectangular 20cm / 30cm [8" / 12"]

* The specifications are subject to change without prior notice.

